

JUN 9 1969

JENT-RE 66-0318

*Reprinted from*

# **ANALYTICAL CHEMISTRY**

---

## **Pulseless High-Pressure Pump for Liquid Chromatography**

**R. E. JENTOFT and T. H. GOUW**

**Chevron Research Co., Richmond, Calif.**

---

**Volume 38, Number 7**

**Pages 949-950, June 1966**

Copyright 1966 by The American Chemical Society and reprinted by permission of the copyright owner.



## Pulseless High-Pressure Pump for Liquid Chromatography

R. E. Jentoft and T. H. Gouw, Chevron Research Co., Richmond, Calif.

**I**N HIGH-PRESSURE liquid chromatography, the eluting solvent is generally pumped into the column by a reciprocating piston or bellows pump. Pulses generated by this pump can often be discerned by the detector. Besides introducing difficulties in the interpretation of the detector output, pulses also have a detrimental effect on the separating efficiency of the column. A major portion of these pulses can be eliminated by the use of surge tanks—e.g., spring-loaded bellows connected to tees in the line between the pump and the column inlet. But even careful design will not remove all pulses from the flow pattern and rather complicated and/or large surge tanks, with concurrent large holdup of material, are necessary to reduce the level of the flow variations to a negligible level. In addition, damage to this pump and surge system is feasible if the outlet line is accidentally blocked.

With the pump described here, high pressures up to 1000 p.s.i. can be generated at the outlet of the pump. In comparison, many of the commercially available pumps are incapable of operating above 100-p.s.i. pressure. Delivery into the chromatographic system can be carried out without any discernible pulses for any length of time and any amount of solvent. Accidental blocking of the outlet line does not damage the pump at all. Construction is simpler than reciprocating plunger pumps of equivalent rating.

### DESCRIPTION AND OPERATION OF PUMP

A schematic drawing of the unit is given in Figure 1. Four 2-inch by 8-inch stainless steel reservoirs, constructed from pipe nipples, are half-filled with mercury, the displacing agent. High pressure air is used as the source of operating pressure. The use of mer-

cury as the displacing agent prevents the dissolution of the air in the solvent.

The outlet of the pump is on the top of Vessel A. For operation up to 250 p.s.i., a dual back-reference flow controller with integral micrometering valves (Millaflow No. 9986 220) is mounted between this outlet and the inlet of the chromatographic column. Because these type flow controllers are currently not available for higher pressure ratings, a fine metering valve is used for work at higher pressures.

The pressure-reducing valve in the air line is adjusted to yield a pressure about one and one-half to two times the desired inlet pressure to the column, but not less than 100 p.s.i. to minimize the effect of the variations in the mercury levels. Prior to introduction into the pump, the solvent in the external reservoir is heated by an infrared lamp. This degases the liquid sufficiently to

prevent air from coming out of solution in the column and creating artifacts in the chromatogram. At the normal mode of operation, neither of the three-way solenoid valves is energized. The air pressure, which is now connected to B, is transmitted through the mercury to the solvent in A. There is vacuum on D, and C is being refilled from the external solvent reservoir. Mercury rises in A as solvent is delivered to the chromatographic column. As soon as the mercury in Vessel A comes in contact with the upper electrode, the solenoids are energized. B is now connected to atmospheric pressure or to a regulated pressure somewhat lower than the operating air pressure in D. The latter is the desired mode of operation at high pressures. The air pressure, which is now on D, refills A from C and, at the same time, keeps up the uninterrupted delivery of the solvent at the




Figure 1. Pulseless high pressure pump





